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**VWR International bvba / spr
B-3001 Leuven**

ENGLISH

You will find the Operating Instructions in English language beginning at page 3.

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Dear Customer

Thank you for purchasing a USC ultrasonic cleaning unit from VWR International. The product you have purchased is one of the most innovative in ultrasonics and provides a maximum of performance and comfort.

This documentation was created in order to make all information concerning the installation, handling and operation of ultrasonic cleaning unit available to the operator. This documentation corresponds to the latest release of the date of publication and is not included in revision service. We reserve the right to make technical changes and improvements at any time.

All devices of the USC-series comply with the regulatory affairs and the normative standards. This is indicated by the **CE**-symbol on the label.

Any enclosed documentation of other vendors serves only for information. We assume no liability for its completeness or correctness.

If you require information extending beyond this documentation, we will be pleased to offer you technical advice at any time.

Technical service

Web Resources

Visit the VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us for information or technical assistance contact your local VWR representative or visit www.vwr.com.

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment. For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you.

VWR International


1. SAFETY INFORMATION.....	5
2. GENERAL.....	8
2.1. Intended purpose / lifetime	8
2.2. Cleaning liquid	8
3. START-UP.....	8
4. OPERATION.....	9
4.1. Operating and display elements	9
4.2. Ultrasonic Cleaning.....	11
4.3. Emptying the tank	11
4.4. Accessories	12
5. MAINTENANCE.....	12
6. APPENDIX.....	13
6.1. Technical Data.....	13
6.2. Operating elements / front view	15
6.3. Connections.....	15
6.4. Dimensions and weights.....	15
6.5. Persistency of the front overlay to solvents	16
6.6. Error messages	17
7. NOTES.....	18

1. SAFETY INFORMATION

Explanation of the used symbols:

Symbol	Explanation
 or 	<p>This symbol points to some danger for persons or the unit, maximum attention is needed in this case. Read the referring chapters especial carefully and follow the requirements strictly.</p> <p>Globally this symbol is also used as reference to follow the operating instructions resp. special advices.</p>
	<p>This symbols shows, that the unit can have hot surfaces.</p>
 or 	<p>Information on the product, that the operating instructions are a part of the unit and must be read in all details.</p>
	<p>CE-sign (Communauté Européenne). The product which carries this symbol meets the requirements of the referring EU-directives (European standards).</p>
	<p>This symbols shows, that the unit must not be disposed with the domestic waste, but collected at special collection points for disposal.</p>
	<p>This symbol marks the unit as unit protection class 1. All electrical conductive parts of the housing are connected to the protective earthing system of the installation, which has earth potential.</p>
	<p>This symbol indicates that the unit has a heating system to warm up the cleaning liquid.</p>

Safety Information

Following safety information must be read and respected before start-up and operation of the ultrasonic cleaning unit:

- Any person using this product must have read and understood this manual. Keep this manual close to the device; it must be permanently accessible to the operator. This will prevent operator error and ensure operation free of trouble.
- Check your device for any transport damage. If you find shipping damage, the device must not be connected to the mains.
- The power cord must be removed from the socket only by pulling the plug, do not pull at the cable. The plug must not be touched with wet / damp hands.
- Power cords must be routed without loops and must not be bent or wrapped around sharp edges.
- The intended use of the ultrasonic cleaning unit is described in this documentation and must be followed.
- Insufficient information about handling or operation can cause irreparable damage and injuries. As a result, the unit must be operated by trained personnel only.
- The ultrasonic cleaning unit must be located on a dry, hard and flat horizontal surface.
- The ultrasonic cleaning unit must be connected to the supply voltage specified on the rating label only. The mains connection must be designed for the rated power consumption.

Always use the correct power cord or –adapter for the available socket outlet.

- Only use cleaning liquids which are approved for use with ultrasonic energy.
- Cleaning liquid in use becomes hot. When cleaning temperature sensitive parts check temperature of cleaning liquid frequently.
- Flammable or explosive liquids and liquids which will corrode the tank material (CrNi steel 1.4301) must not be used. Liquids which contain halogens (chlorine, bromine or iodine ions) can cause irreparable damages.
- Never switch on the ultrasonic cleaning unit without cleaning liquid. The cleaning tank must be filled at least to 7/8 of volume (approx. up to the seam). Also ensure that the tank is filled enough during operation.
- The temperature rise of the cleaning liquid caused by heaters or ultrasonic energy can result in steaming. Ensure adequate ventilation at the environment of installation.
- Fluid must not get inside the unit. Wipe off any fluid which has slopped over and is leaking down the side walls as soon as possible
- The items to be cleaned must not be placed directly on the bottom of the tank. Use a suitable cleaning basket.
- Never leave the ultrasonic cleaning unit unattended during operation.
- If the unit does not have to be ready for operation, disconnect it from the power supply (unplug it from mains).
- Never operate the unit if it has been damaged. Inspect the device and the power cord for outer damages frequently. Repairs must only be carried out by the manufacturer.
- Observe the statutory regulations when disposing of the cleaning liquid and disposing the unit.

WARNING:

To avoid risk of **electric shock**, this equipment must only be connected to supply mains with **protective earth**.

Take special note on using the **proper power cable and/or adapter** for the available power outlet.

WARNING:

There is a danger of **injury and burns** through:

- Ultrasonic energy!
Never grab in the tank with your hands when the ultrasonic energy is switched on. Always switch off the ultrasound when inserting or removing the items to be cleaned!
- Hot liquid!
- Hot surfaces!
- Toxic vapor!

During operation it is permitted to touch or to handle the unit with the handles or at the control panel only.

After use the unit must cool down min. 30min prior to touching the housing.

WARNING:

For usage of the cleaning unit and handling with cleaning agent you generally have to wear adequate personal protective equipment (PPE) corresponding to the standards:

- Safety gloves
- Safety glasses
- Ear protection
- Breathing protection
- Protective clothing
- Protective boots

WARNING:

No modification of this equipment is **allowed**.

Do not carry out any repair. In case of service, send the unit back to the manufacturer.

General / Start-up

2. GENERAL

2.1. Intended purpose / lifetime

The USC ultrasonic cleaning units were especially developed for cleaning tasks in laboratory facilities of commercial use.

The expected life time of the USC ultrasonic cleaning units depends basically on the quality of maintenance and cleaning. With proper use, maintenance, cleaning and any necessary repairs a lifetime of five years can be expected. This does not include consumables such as cable / connectors, cleaning baskets and similar.

Only accessories which are approved by or can be obtained from the manufacturer have to be used.

2.2. Cleaning liquid

The cleaning liquid which has to be used depends on the material of the items to be cleaned and the type of soiling. Please contact the supplier of the cleaning agent.

3. START-UP

- Check if the mains connection, which must have a protective earth conductor, confirms with the specification on the rating label of your device.
- Connect the power socket (circular connector at the cable) to the flange connector at the rear of the unit and lock the connection by turning the cap nut. The cap nut must latch noticeably and needs to be disconnected in locked condition with increased force only.

If the power cord is already installed to the unit, check if it is locked properly as described before.

- The ultrasonic cleaning unit must be located on a dry, hard and flat horizontal surface.
- Fill the cleaning tank with a suitable cleaning liquid depending on the application. The level must be at least 7/8 of volume (approx. up to the seam).

Take note of the respective safety regulations for handling chemicals.

- Connect the ultrasonic cleaning unit to the mains, when the display lights up, the unit is ready for operation.

WARNING:

To avoid risk of **electric shock**, this equipment must only be connected to supply mains with **protective earth**.

Take special note on using the **proper power cable and/or adapter** for the available power outlet. As a result of the design, the EU connector (CEE 7/7) fits mechanically into the DK socket outlet (Denmark) but does not connect the protective earth conductor (PE).

**Never use the EU connector in combination with DK socket outlets
DANGER TO LIFE AND LIMB!**

Adapters are one-time mountable only and cannot be removed after applying to the connector. Therefore be sure to mount the proper adapter.

WARNING:

The power plug is used to disconnect the unit from mains.

The used **power outlet** must be **permanently accessible** and must not be covered by the unit or any other objects, so the mains plug can be unplugged at any time.

4. OPERATION

The units have a digital operator interface which allows comfortable setting and operation.

Depending on the version there are different configurations available:

Configuration	Series			
	T	TH	THD	THD/HF
Table top cleaning unit with digital control and display, ultrasonic frequency approx. 45kHz	X	X	X	
Table top cleaning unit with digital control and display, ultrasonic frequency approx. 132kHz				X
Timer for ultrasonic operation adjustable from 1 to 99 minutes in steps of 1 minute or continuous	X	X	X	X
Integrated heater, liquid temperature adjustable in steps of 1K up to 80degC		X	X	X
Ultrasonic power adjustable in nine steps from approx. 40% to 100 % power			X	X
DEGAS-function			X	X

DEGAS-function:

For degasification of the cleaning fluid the unit can be operated in a so-called DEGAS mode. The cleaning liquid will be also degassed during standard cleaning operation but this can be considerably speeded up in the DEGAS mode. With degassing the cleaning result during the following cleaning process is considerably more effective and uniform.

During the DEGAS mode, the ultrasound is automatically switched on and off in short intervals. This causes the entrapped air to form into larger bubbles during the acoustic irradiation of the fluid. The bubbles can then rise freely to the surface in the pauses between the acoustic irradiation.

4.1. Operating and display elements

Button 'Sonic / Ultraschall / Ultrasons'

With this button you can switch on and off the ultrasonic operation. The operating status 'Ultrasonic ON' is indicated by the green light-emitting diode (LED) in this button.

Button 'Heater / Heizung / Chauffage' (TH, THD and THD/HF types only)

With this button you can switch on and off the liquid heater. The operating status 'Heating ON' is indicated by the green LED in this button.

Button 'Degas / Dégazage' (THD and THD/HF types only)

With this button you can switch on and off the DEGAS mode (see description below). The operating status is indicated by the green LED in this button.

Button Taste 'Select / Auswahl / Choix'

With this button the following operating parameters can be selected and set (time, temperature, power):

Operation

Display 'Time / Zeit / Durée' (two digits)

During ultrasonic **switched off**, the pre-selected cleaning time is shown in minutes in this display. When the ultrasonic is **switched on**, the remaining cleaning time is shown and the decimal point in the display flashes, thus indicating that the pre-selected cleaning time is counting down.

Display 'Temperature / Temperatur / Température' (two digit, TH, THD and THD/HF types only)

In cleaning mode, the current bath temperature is shown in degC in this display. If the pre-selected temperature has not yet been reached, the decimal point in the display flashes; the bath is being heated. Once the pre-selected temperature has been reached, the decimal point is steady. During the process 'Setting the bath temperature', the target temperature for the cleaning bath is shown in the display.

Display 'Power / Leistung / Puissance' (single digit, THD and THD/HF types only)

The current set value of the ultrasonic power is shown in steps from 1-9 in this display. Step 9 corresponds to maximum power (100%).

Setting the cleaning time

With ultrasonic switched off press the 'Select / Auswahl / Choix' button until the 'Time / Zeit / Durée' display flashes. The desired cleaning time can be set in steps from 1 – 99 minutes with the '+' and '-' adjusting buttons.

The time pre-selection can be disabled by pressing the '-' button until two horizontal lines are shown in the display ('- -').

Once the desired time has been set, release the '+' and '-' adjusting buttons and wait until the display stops flashing (approx. 5 sec). The value is stored now.

Setting the bath temperature ('TH', 'THD' und 'THD/HF' types only)

With ultrasonic switched off press the 'Select / Auswahl / Choix' button until the 'Temperature / Temperatur / Température' display flashes. The desired liquid temperature can be set up to 80degC with the '+' and '-' adjusting buttons.

Once the desired temperature has been set, release the '+' and '-' adjusting buttons and wait until the display stops flashing (approx. 5 sec). The value is stored now.

Setting the ultrasonic power ('THD' und 'THD/HF' types only)

The ultrasonic power can be set during ultrasonic on or off.

Ultrasonic 'OFF':

During ultrasonic switched off press the 'Select / Auswahl / Choix' button until the 'Power / Leistung / Puissance' display flashes. The desired power can be set in steps from 1 – 9 with the '+' and '-' adjusting buttons (step 9 = max. power).

Once the desired power is set, release the '+' and '-' adjusting buttons and wait until the display stops flashing (approx. 5 sec). The value is stored now.

Ultrasonic 'ON':

During ultrasonic switched on press the 'Select / Auswahl / Choix' button until the 'Power / Leistung / Puissance' display flashes. The desired power can be set in steps from 1 – 9 with the '+' and '-' adjusting buttons (step 9 = max. power).

Once the desired power is set, release the '+' and '-' adjusting buttons and wait until the display stops flashing (approx. 5 sec). The value is stored now.

Degassing the fluid (DEGAS, 'THD' und 'D/HF' types only)

In order to obtain an optimum cleaning process after the unit has been refilled, first of all the liquid should be degassed. The cleaning result is considerably impaired by the solved gas bubbles (usually ambient air) in the liquid. Ultrasonic vibrations dispense these gases.

In the DEGAS mode the unit is periodically switched on and off, the gas bubbles can raise out of the fluid in the OFF period.

Degassing can be speeded up by adding a few drops of a surface-active agent (e.g. washing-up liquid) to the cleaning fluid.

Most ultrasonic cleaning agents already contain a surface-active agent.

The degassing process does not have to be repeated when reusing the cleaning liquid.

Heating the cleaning bath (TH, 'THD' und 'THD/HF' types only)

If a specified bath temperature is required for the cleaning operation, the cleaning bath can be heated up to a pre-selected temperature (see process 'Setting the bath temperature'). Cleaning is normally considerably more effective if the cleaning fluid is warmed up.

WARNING:

Consider the **temperature tolerance** of the material to be cleaned!

4.2. Ultrasonic Cleaning

- Immerse the items to be cleaned in the cleaning bath the distance of the items from the bottom must be at least 20 mm (0.8 inches).
If required, use the cleaning basket which is available as an accessory.

- Set the required cleaning time, temperature and ultrasonic power for the cleaning process (if applicable).

The parameters 'required cleaning time' and 'ultrasonic power level' excessively depend on the cleaning fluid, temperature, type and degree of soiling of the items requiring cleaning. The time can vary from a few seconds to a few minutes.

The best combination of the parameters may have to be experimentally determined.

- Start the cleaning process by switching on the ultrasonic energy.
- After the cleaning process remove the items out of the liquid with suitable care. We recommend rinsing the parts in clean water.
- Disconnect the ultrasonic cleaning unit from mains if the unit is not needed.

WARNING:

Please note that the bath temperature rises considerably with duration of the ultrasonic operation. For safety reasons the heater and the ultrasonic energy are switched off if the bath temperature exceeds approx. 90degC (194F), an error message appears in the display.

Never put the items to be cleaned directly on the bottom of the tank.
Never clean too many items at once, the cleaning result can deteriorate.

4.3. Emptying the tank

If the cleaning result becomes worse or if the cleaning liquid has got contaminated, the cleaning liquid should be replaced.

Prior to any emptying or cleaning the unit must be disconnected from mains.

For this purpose the units USC 900, USC 1200, USC 1700, USC 2100 and USC 2600 have a drain on the rear side of the unit.

4.4. Accessories

- **Tank cover (stainless steel)**
To prevent evaporation of the tank contents and odor nuisance.
- **Cleaning basket (stainless steel)**
For carrying the items to be cleaned. The basket can be suspended in the cleaning tank and it also can be placed on the edge of the tank with the handles for draining.
- **Cover with openings for beakers (different sizes)**
 - easy cleaning of small parts
 - less needs of cleaning solvent
 - cleaning in different cleaning solvents

The ultrasonic energy is coupled inside the filled beaker by the liquid around through the glass. Fill the tank with water and add a surface-active agent (e.g. washing-up liquid).

5. MAINTENANCE

WARNING:

The ultrasonic cleaning unit needs low maintenance.
Regard following notes:

- Prior to any emptying or cleaning the unit must be disconnected from mains.
- Keep the device clean. Wipe the housing frequently with a damp cloth. Remove cleaning agent residues immediately.
- Ensure that the cleaning liquid does not soil the housing during cleaning.
Clean, replace and dispose the cleaning liquid in accordance with the supplier's recommendations.
- Avoid scratching the inner tank surface, this accelerates cavitation (erosion). Ensure the items to be cleaned have no direct contact with the bottom or inner sidewalls of the tank and try to keep none or only small amounts of pollution on the bottom of the cleaning tank. If necessary, the cleaning liquid should be drained off prematurely, so the tank can be cleaned.
-  Check the unit and the power cord for damages frequently and prior to any use.
A visibly damaged unit or cable must not be put in operation and must be disconnected from the mains.
Do not carry out repairs yourself. Send the unit to the manufacturer for servicing or repair.
- There are any parts to maintain or which have to be changed by the user inside the unit.

If the power cord is damaged, it can be removed by unlocking the cap nut (circular plug at rear side of unit). The power cord must be replaced with the manufacturer's original part only.

6. APPENDIX

6.1. Technical Data

Intended Use Cleaning applications in laboratory facilities of commercial use.

Executions Table top cleaning unit with digital control and display

Unit series	Execution
USC - T	<ul style="list-style-type: none"> ultrasonic energy with approx. 45kHz
USC - TH	<ul style="list-style-type: none"> ultrasonic energy with approx. 45kHz built-in liquid heater
USC - THD	<ul style="list-style-type: none"> ultrasonic energy with approx. 45kHz built-in liquid heater adjustable ultrasonic power DEGAS-function
USC - THD/HF	<ul style="list-style-type: none"> ultrasonic energy with approx. 132kHz built-in liquid heater adjustable ultrasonic power DEGAS-function

Output ultrasonic energy with frequency of approx. 45kHz resp. 132kHz
switch-on duration: 100% (continuous operation)

Environment conditions

Operating ambient temperature: 0degC ... +40degC (32F ... 104F)
 Transportation & storage temperature: - 20degC ... +65degC (- 4F ... 149F)
 Humidity: 20% ... 70%, non-condensing
 Altitude (sea level): max. 2000m (a.s.l.)

Classification:

Protection class: class I equipment (earthed housing)
 Protection category according IEC 529: IP 21
 Moisture degree 2
 Over voltage category II

Mains connection 230V-versions:

AC230V ±10%, 47...53Hz, 1~, N, PE / 16A

Unit type (230V)		HF-power		connected load		heater power *)
		peak	effective	'T'	'TH', 'THD'	'TH', 'THD'
USC 300	S TH THD	160W	80W	150VA	370VA	200W
USC 500		200W	100W	170VA	390VA	200W
USC 600		240W	120W	190VA	620VA	400W
USC 900		400W	200W	320VA	950VA	600W
USC 1200		360W	180W	280VA	910VA	600W
USC 1700		480W	240W	360VA	1170VA	800W
USC 2100		480W	240W	360VA	1170VA	800W
USC 2600		600W	300W	420VA	1430VA	1000W

Unit type (230V)		HF-power		connected load	heater power *)
		peak	effective		
USC 300	THD/HF	160W	80W	400VA	200W
USC 600		240W	120W	660VA	400W
USC 1200		320W	160W	960VA	600W
USC 2100		420W	210W	1250VA	800W

*) approx. values

Mains connection

AC 115V ± 10%, 47...63Hz, 1~, N, PE / 15A

Unit type (115V)		HF-power		connected load		heater power *)
		peak	effective	'T'	'TH', 'THD'	'TH', 'THD'
USC 300	S TH THD	160W	80W	160VA	370VA	200W
USC 500		200W	100W	180VA	390VA	200W
USC 600		240W	120W	190VA	620VA	400W
USC 900		400W	200W	320VA	950VA	600W
USC 1200		360W	180W	280VA	910VA	600W
USC 1700		480W	240W	370VA	1180VA	800W
USC 2100		480W	240W	370VA	1180VA	800W
USC 2600		600W	300W	430VA	1440VA	1000W

Unit type (115V)		HF-power		connected load	heater power *)
		peak	effective		
USC 1200	THD/HF	320W	160W	980VA	600W
USC 2100		420W	210W	1280VA	800W

*) approx. values

6.2. Operating elements / front view


Series 'T'


Series 'TH'

Series 'THD'

Series 'THD/HF'

6.3. Connections

Mains connection, single phase 230V:

- power cord detachable with power plug CEE 7/7
- cord length approx. 2.0m (6.5ft)

Adapters available for

- United Kingdom – BS-1363
- Switzerland – SEV1011
- Denmark – DS 60884-2-D1
- China CPCS-CCC / Australia – AS 3112

Attention:

Adapters are one-time mountable only and cannot be removed after applying to the connector. Therefore be sure to mount the proper adapter.

Mains connection, single phase 115V:

- power cord detachable with power plug NEMA 5-15P
- cord length approx. 2.0m (6.5ft)

6.4. Dimensions and weights

Unit type	Tank capacity (liters)	inner dimensions (LxWxD / mm)	outer dimensions (LxWxD / mm)	weight (kg)
USC 300	2.8	237x134x100	265x162x235	4.1
USC 500	4.2	237x134x150	265x162x295	4.9
USC 600	5.4	297x148x150	325x176x295	5.5
USC 900	9.2	502x134x150	515x150x270	8.5
USC 1200	12.3	297x237x200	325x265x335	8.2
USC 1700	16.8	324x297x200	352x325x335	9.7
USC 2100	19.6	502x297x150	530x325x365	12.7
USC 2600	26.1	502x297x200	530x325x365	12.9

Appendix

6.5. Persistency of the front overlay to solvents

The front overlay is based on a polyester foil with a biaxial alignment and has therefore a higher persistence against solvents compared with e.g. polycarbonate or PVC. It is stronger and more durable than other foils generally used for plastic foil keyboards and front panels.

The material is resistant against following chemicals over a period of more than 24 hours without any perceptible changes:

Ethanol Cyclohexane Diaceton alcohol Glycol Isopropanol Glycerin Methanol Triacetine Dowandol DRM/PM	Formaldehyde 37 – 42 % Acetaldehyde Aliphatic Toluene Xylene Thinning agent (white spirit)	Trichloroethylene Ethyl acetate Diethyl ether N-butyl acetate Amyl acetate Butylcellosolve Ether
Acetone 2-butanone Dioxane Cyclohexane MIBK Isophorone	Formic acid < 50 % Acetic acid < 50 % Phosphoric acid < 30 % Hydrochloric acid < 36 % Nitric acid < 10 % Trichloroacetic acid < 50 % Sulphuric acid < 10 %	Soda chlorine < 20 % Hydrogen peroxide < 25 % Potash soft soap Detergents Tensides Softening agents Iron (II) chloride (FeCl ₂) Iron (III) chloride (FeCl ₃) Dibutyl phthalate Dioctyl phthalate Sodium carbonate
Ammonia < 40 % Sodium hydroxide solution < 40% Potassium hydroxide Alkali carbonate Dichromate Potassium hexacyanoferrate Acetonitrile Sodium bisulphate	Cutting emulsions Diesel oil Varnish Paraffin oil Castor oil Silicone oil Turpentine substitute Brake fluid Decon Aviation fuel Petrol Water Salt water	

E

6.6. Error messages

The ultrasonic cleaning units have an error detection system which shows the following messages in the two most left displays in case of an error:

- **nC or E0**
In this case, an internal error occurred.
Disconnect the ultrasonic cleaning unit from mains by unplugging the power connector and reconnect it. Restart the cleaning operation, if the error message still occurs, the device must be sent to the manufacturer for servicing.
- **E1**
In this case, an error has been detected in the ultrasonic generator or the generator control system.
Disconnect the ultrasonic cleaning unit from mains by unplugging the power connector and reconnect it. Restart the cleaning operation, if the error message still occurs, the device must be sent to the manufacturer for servicing.
- **E2**
In this case, the maximum permissible bath temperature of 90degC has exceeded. The heater and the ultrasonic operation are automatically switched off.
Press the button 'Select / Auswahl / Choix' to quit the error condition and wait minimum 5 minutes before restarting the cleaning operation.
- **E3**
In this case, the temperature protection of the ultrasonic generator has automatically switched off the unit.
Please check that there is enough liquid in the cleaning tank (7/8 of the volume), the ambient temperature is in the specified range and there is sufficient ventilation of the unit. Press the button 'Select / Auswahl / Choix' to quit the error condition and wait minimum 5 minutes before restarting the cleaning operation.

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VWR Singapore Pte Ltd
VWR International Oy
VWR International AB
VWR International Ltd