

Faculty / Service Area:	Faculty of Health Sciences and Sport	Location:	Cottrell 2B145						
Description of work task / equipment /a	rea being assessed:								
Body composition analysis using a Dual e	energy x-ray absorptiometry scanner								
Head of division	Prof Jayne Donaldson	Safety office	r	Dr Nidia Rodriguez Sanchez					
Completed by:	Frank Kelly	Date: (01/10/2010					
Reviewed by (Line Manager):	Chris Grigson	Date:		10/07/2023					
		Date of next review:		01/08/2024					
Equipment used	GE Medical Lunar iDXA Scanner	GE Medical Lunar iDXA Scanner							
Categories of people involved	Staff, UG, PG, Visitors								
Duration of activity	30 mins	Frequency of activity		Patients < 10 scans per year Individual staff conduct < 150 scans per year					
Legal compliance to standards and regulations required	The Ionising Radiations Regulations 2017 (IRR) <u>https://www.legislation.gov.uk/uksi/2017/1075/contents/made</u> Health and Safety at Work act 1974 (HASAWA) <u>https://www.hse.gov.uk/legislation/hswa.htm</u> Management of Health and Safety at Work Regulations 1999 (MHSWR) <u>https://www.legislation.gov.uk/uksi/1999/3242/contents/made</u>								

Risk assessment RA10

https://sportsciencesafety.stir.ac.uk



			Provision of Work Equipment Regulations 1998 (PUWER) https://www.hse.gov.uk/work-equipment-machinery/puwer.htm The Control of Substances Hazardous to Health Regulations 2004 (COSHH) https://www.hse.gov.uk/coshh Personal Protective Equipment at Work (Amendment) Regulations 2022 https://www.legislation.gov.uk/uksi/2022/8/contents/made Data Protection Act 2018 https://www.legislation.gov.uk/ukpga/2018/12/contents/enacted Access to Health Records Act 1990 https://www.legislation.gov.uk/ukpga/1990/23/contents							
Change log			25/10/2 03/08/2 05/12/2 13/12/2 10/07/2	2016 2022 Change of Safe Change of Mer New format 2022 V1.0 Improve gram Add reference 2022 V1.1 Added data pr Listed associat Subjects are p 2023 V1.2 Added hyperli	ety officer dical praction s otection ced docume atients nks	tioner ents				
What are the Hazard Who might be hazards? category harmed and how?		t be how?	What are you already doing to control the risks?	Risk rating	What additional controls (if any) are required to reduce the risks?	Risk rating	Action by who?	Action by when?	Date of completion	
DXA scanning (X-rays)	F4	All Exposure to ic radiation	onising	Radiation risk assessment: 221205_DXARadiationRA.docx	Medium	Operator training regularly refreshed				



	DVA installed in dedicated room				
	with suitable warning signs and				
	access control.				
	Local Rules displayed in				
	supervised area:				
	220802 localDLUESCottroll2D14E				
	All staff operating in the				
	laboratory IR(MR)R Trained and				
	Trained by the equipment				
	manufacturer				
	Standard operating procedure				
	Monthly inspection				
	Monthly hispection				
	A second s				
	Annual maintenance by				
	manufacturer				
	Operators have procedures to				
	protect them from the dangers				
	of ionising radiation (safe				
	distance from DXA scanner				
	marked out). No personal				
	protoctive equipment				
	required.				
	Passing foot traffic controlled by				
	warning signs placed on door of				
	DXA room to prevent accidental				
	exposure. Warning signs light to				
	indicate a scan is in progress				
		1			



		(At 1 m dose rate is <10 μSv/h – information from supplier ^a)				
	Patients Uncontrolled exposure to ionising radiation	Referral process used to assess the patient's requirements for DXA scanning, if patient is not suitable or needs do not outweigh the risks then they are not scanned. Decision to screen lies with the medical practitioner (Dr Chris Kelly)	Medium			
		Patients identity established beyond doubt by operator (see procedures) so as to ensure that the wrong person is never scanned.				
		Informed consent forms used to inform Patients of risks. Withdrawal from scan if patient doesn't consent.				
	Female Patients Exposure to ionising radiation of pregnant women.	Pregnancy is screened for, as exposure to ionising radiation can cause harm in the case of pregnancy existing. Patients are asked about possibility of them being pregnant, if possibility of	Medium			



	pregnancy exists then they will				
	not be tested (until proven other				
	wise)				
	,				
	All Patients of child bearing age				
	to submit a urine sample to				
	establish a negative pregnancy				
	test before scan.				
Patients	Operators work within the	Modium			
	procedures to ensure scans are	Medium			
Unnecessary	carried out correctly and do not				
exposure to ionising	need to be repeated, thus				
radiation	Patients do not undergo				
	unnecessary scans.				
	Quality control procedures are in				
	place to ensure equipment is				
	working properly and scans do				
	not need to be repeating, thus				
	Patients do not undergo				
	unnecessary scans.				
All	Operators have procedures to	High			
	protect them from the dangers	Ingri			
Overexposure to	of ionising radiation (safe				
ionising radiation	distance from DXA scanner is				
(Routine use)	marked out). No personal				
	protective equipment required.				
	Radiation safety assessment				
	carried out to determine safety				
	of equipment and safe operating				
	distances for operators and				



room.							
DXA fails to de-energise at end of scan. DXA serviced in accordance with manufacturer's recommendations.	High	Contingency plan in Local Rules					
Phantom stored in special case, care taken by operators when moving it	Low						
Scanner arm should be moved via provided controls to suitable position to allow patient access. See SOP	Low						
Scanner arm moves over patient head, if patient were to sit up risk of injury may occur. Operators instruct all patients to remain still throughout scan and explain all the procedures to all patients prior to scans	Low						
	DXA fails to de-energise at end of scan.DXA serviced in accordance with manufacturer's recommendations.Phantom stored in special case, care taken by operators when moving itScanner arm should be moved via provided controls to suitable position to allow patient access. See SOPScanner arm moves over patient head, if patient were to sit up risk of injury may occur. Operators instruct all patients to remain still throughout scan and explain all the procedures to all patients prior to scans commencing See SOP	DXA fails to de-energise at end of scan.HighDXA serviced in accordance with manufacturer's recommendations.HighPhantom stored in special case, 	DXA fails to de-energise at end of scan.HighContingency plan in Local RulesDXA serviced in accordance with manufacturer's recommendations.HighContingency plan in Local RulesPhantom stored in accordance with moving itLowImage: Contingency plan in Local RulesScanner arm should be moved via provided controls to suitable position to allow patient access. See SOPLowScanner arm moves over patient head, if patient were to sit up risk of injury may occur. Operators instruct all patients to remain still throughout scan and explain all the procedures to all patients prior to scans commencing. See SOPLow	DXA fails to de-energise at end of scan.HighContingency plan in Local RulesDXA serviced in accordance with manufacturer's recommendations.HighContingency plan in Local RulesPhantom stored in accordance with moving itLowLowScanner arm should be moved via provided controls to suitable position to allow patient access. See SOPLowScanner arm moves over patient head, if patient were to sit up risk of injury may occur. Operators instruct all patients to remain still throughout scan and explain all the procedures to all patients prior to scans rommering. See SOPLow	DXA fails to de-energise at end of scan.HighContingency plan in Local RulesDXA serviced in accordance with manufacturer's recommendations.HighContingency plan in Local RulesPhantom stored in special case, care taken by operators when moving itLowImage: Construct of the second seco	DXA fails to de-energise at end of scan. High Contingency plan in Local Rules DXA serviced in accordance with manufacturer's recommendations. High Contingency plan in Local Rules Phantom stored in special case, care taken by operators when moving it Low Low Scanner arm should be moved via provided controls to suitable position to allow patient access. See SOP Low Low Scanner arm moves over patient head, if patient were to sit up risk of injury may occur. Low Low Scanner arm moves over patient head, if patient were to sit up risk of injury may occur. Low Low Scanner arm moves over patient head, if patient were to sit up risk of injury may occur. Low Low Scanner arm moves over patient head, if patient were to sit up risk of injury may occur. Low Low Scanner arm moves over patient head, if patient stor remain still throughout scan and explain all the procedures to all patients prior to scans commercing. See SOP Low	DXA fails to de-energise at end of scan.High HighContingency plan in Local RulesDXA serviced in accordance with manufacturer's recommendations.HighContingency plan in Local RulesPhantom stored in special case, care taken by operators when moving itLowImage: Control of the second seco

Physiology, Exercise & Nutrition Research Group



Unintended release or loss of patient health data ^c	F4	Patients Patients can no longer access their data as required by the DPA and AHRA	Operators store only anonymised results on secure platforms. All patient identification lists are kept separately on secure platforms.	Low			
		Patient data is made publically available leading to harm The University of	The safety officer is responsible for the safe archiving of historic data and patient identification lists. All researchers must securely pass their data to the				
		Stirling The University is shown to have broken data protection law which results in prosecution and loss of reputation	safety officer at the end of their study Scan data is backed up to an encrypted hard drive once a month				
Deferences		The University is subject to litigation through the civil courts					

a Lunar encore, Safety and Technical Specification Manual, GE Healthcare, rev 3, May 2009 https://www.gehealthcare.com/-/media/20fc07d1369e4d15acae5732090559db.pdf

b Work with ionising radiation, Ionising Radiations Regulations 2017, Approved Code of Practice and guidance, Jan 2018 https://www.hse.gov.uk/pubns/priced/l121.pdf



c Information Commisisioner's Office, Guide to the UK General Data Protection Regulation (UK GDPR), https://ico.org.uk/for-organisations/guide-to-data-protection/guide-to-the-general-data-protection-regulation-gdpr